Vertex Castaquick

EN For all types of removable denture repairs

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EN Instruction For Use

Introduction

The following instructions for use are for dental professionals who use Vertex Castaquick to repair removable full and partial denture bases. This material is intended exclusively for professional dental work. This instruction for use also provides information about safety and environmental aspects. In case more information is needed contact the reseller.

Intended use

Vertex Castaquick is used for the repair and/or relining of all types of removable full and partial denture bases. The product is CE certified as a Class IIa medical device. The technical instructions provided in this user guide should be followed. Any deviations from these instructions may have a negative effect on the intended result and will not guarantee the quality of the result.

Description & Effects

Vertex Castaquick is a 2-component acrylic system. The acrylic system is formed by a powder (polymer) and a liquid (monomer). Vertex Castaquick powder is used with Vertex Castavaria liquid. The combination of powder and liquid is converted into a hard finished product by a self-initiating process that is accelerated by heat.

Contra-indications

Vertex Castaquick should not be used for any other purpose than for the repair and /or relining of removable full and partial denture bases. Any deviation from this instruction for use may have an adverse effect on the chemical and physical quality of Vertex Castaquick. In case of an allergic reaction, please contact a medical physician.

Vertex Castaquick Powder component Hazard & Precaution

Warning

Hazard statements
H411 - Toxic to aquatic life with long lasting effects.

Precautionary statements

P273 - Avoid release to the environment.
P391 - Collect spillage.
P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

EUH Statements

EUH208 - Contains dibenzoyl peroxide; benzoyl peroxide, methyl methacrylate; methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate. May produce an allergic reaction.

For professional users only.

Vertex Castavaria Liquid component Hazard & Precaution

Danger

Hazardous ingredients

methyl methacrylate; methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate, 2-Propenoic acid, 2-methyl-, 1,2-ethanediy l ester, 2-hydroxyethyl methacrylate, 2-(2H-benzotriazol-2-yl)-p-cresol

Hazard statements

H225 - Highly flammable liquid and vapour. H315 - Causes skin irritation. H317 - May cause an allergic skin reaction. H335 - May cause respiratory irritation.

Precautionary statements

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 - Wear protective gloves.
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.
P370+P378 - In case of fire: Use foam, dry extinguishing powder, carbon dioxide (CO2) to extinguish.
P403+P233 - Store in a well-ventilated place. Keep container tightly closed.

For professional users only.

Processing

To repair a denture roughen the surface in the immediate vicinity of the expansion or fracture. Make sure the denture base surface is free of dust and grease. To repair denture teeth roughen the base side of

artificial acrylic teeth and treat additionally with the monomer or a contact fluid (e.g. Vertex Acrybond).

The mixing ratio of Vertex Castaquick is 1 ml liquid : 1,7 gr powder. Measure both components in proportion. Pour the liquid into a chemical-resistant bowl first and the powder second. Mix gently for 20 seconds. Cover the bowl and let the mixture rest for 5 - 8 minutes. The mixture will have reached the dough stage at the end of the rest period. The mentioned mixing, dough and working times apply for a material and ambient temperature of 22°C. A lower or higher temperature will have a delaying or accelerating effect on the treating times.

Take the dough out of the mixing bowl and apply it to repair the denture base. Polymerize in a pressure vessel. Polymerize at 55°C at 2½ bar for 10 minutes.

Vertex Castaquick		
Mixing ratio	Mix	Temp
1 ml liquid: 1,7 g powder	20 sec	22°C
Doughtime		Polymerize
max. 8 min		55°C /2½ bar/ 10 minutes

Finishing

Finish and polish the polymerized parts using conventional dental methods and instruments.

Storage conditions, expiry date and transport

Store the polymer in a cool, dry environment. Close the packaging properly after each use. Store the liquid only in the original container in a cool well ventilated place. Keep container tightly closed. Store locked up. Keep in fireproof place. The plastic components have a use-by date indicated on the product label. After the use-by date, the plastic components are no longer guaranteed in terms of treatment. Transportation is restricted by regulations on transporting hazardous materials of the Vertex Castavaria liquid component (UN 1247), Vertex Castaquick powder component (UN 3077).

Waste disposal

Treat remnants and empty liquid packaging as chemical waste. Empty containers of the liquid component should be handled with care because residual vapours are flammable. Special disposal requirements are applicable. Check by your local, federal, or other regulatory agencies for disposal requirements.

Information for dental professionals and patients

The following information should be transferred from dental professional to patient:

- Periodic check-ups by a dentists or dentist are required to monitor changes in the form and structure of the residual ridges of the patient that uses a denture.
- In case of breakage of a denture base, potentially damage the mucosa, the palate or other parts of the patients mouth or esophagus can occur. Contact a dentist.
- In case of an allergic reaction, contact a medical physician.
- Vertex Castaquick repairs and/or relinings can be cleaned with non-aggressive and non-abrasive dental cleaning products.

If you choose not to follow this Instruction for Use, Vertex Dental cannot be held accountable for adverse effects on the safety and performance of the denture base repaired with Vertex Castaquick. Should any serious incident occur, report this to the manufacturer and competent authority.

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