

# Vertex™ Alma Gauge



## Features using the Vertex™ Alma Gauge

- More accurate work due to the exact measurements in the beginning
- Occlusal rims can be trimmed more precise
- The arch can easily be copied
- The setup of the denture can easily be checked with the available readings of the Vertex™ Alma Gauge
- The finished denture can easily be checked with the available readings of the Vertex™ Alma Gauge
- Very convenient in a Class II or Class III situation for the correct relationship
- Fewer remakes which saves time for dentists, dental technicians and patients

**Vertex** Dental



1 Place an accompanied transparent plastic sheet on the Vertex™ Alma Gauge, then place the existing denture on the Vertex™ Alma Gauge. Write down the horizontal and vertical dimension.



2 Use the accompanied fine liner to copy the arch of the denture on the transparent plastic sheet. Send the sheet including the measurements and the impressions to the dental laboratory so they can make a baseplate including a wax rim accordingly.



3 The lab will make a base plate and wax rim according to the received measurements. They can use the Vertex™ Alma Gauge to check the measurements on the baseplate. When the baseplate and the wax rim is made and checked, it can be sent back to the dentist for the final adjustments.



4 When everything is returned with the final individual information from the dentist, the dental technician can start with the setup. When the setup is complete he can check the setup with the Vertex™ Alma Gauge and the obtained measurements.



5 The obtained plastic sleeve is used for the correct form of the arch. Now the denture can be sent to the dentist for the trial fit when everything is approved by the dentist and patient the denture can be sent back to the lab and can be finished by the lab.



6 When the denture is finished there will be a final check for the horizontal and vertical dimensions. When everything is correct the denture can be fitted in the patients mouth.

#### Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist  
P.O. box 10 | 3700 AA Zeist | The Netherlands  
Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188  
info@vertex-dental.com | www.vertex-dental.com

#### Vertex-Dental Asia Pte Ltd

The Alpha #04-16 | Singapore Science Park II  
10 Science Park Road | Singapore 117684  
Tel. +65 6872 9315 | Fax +65 6775 7963  
info.sg@vertex-dental.com | www.vertex-dental.com