

Quint Teeth Mould Chart



Quint Synthetic Resin Teeth - Upper Anteriors

Square



M1
41.3
8.4
8.0



M2
43.5
9.8
8.3



M3
46.0
10.3
8.7



M4
46.3
10.8
8.8

Ovoid



J1
41.8
9.0
7.6



J2
42.7
9.4
7.8



J3
46.3
9.6
8.6



J4
46.4
10.3
8.8



J5
52.0
12.5
9.8

Tapering



R1
41.8
9.9
8.1



R2
43.2
10.5
8.4




R3
47.0
10.7
8.4



R4
48.5
10.9
8.9



R5
50.1
13.0
9.5



All dimensions are in millimeters

Quint Synthetic Resin Teeth - Lower Anteriors



AL1
32.1
8.2
4.7



AL2
34.8
8.5
5.2



AL3
34.6
9.4
5.0



AL4
35.8
9.4
5.7



AL5
37.4
9.0
5.7



AL6
38.1
9.6
5.7



AL7
38.8
11.0
5.6



AL8
39.0
10.7
6.0

Quint Synthetic Resin Teeth - Posteriors

Vertex Dental



07

29.9
8.6
31.8
8.2



08

29.1
8.5
31.6
8.0



09

31.8
9.1
34.3
8.7



10

34.1
9.7
36.6
9.1



11

30.3
8.5
33.5
7.4



12

28.0
7.6
30.0
6.9



13

30.0
9.0
31.5
7.5



14


30.5
9.2
31.2
8.8



15


32.8
10.2
33.5
9.4

08



All dimensions are in millimeters

11



All dimensions are in millimeters

Quint Synthetic Resin Teeth



Combinated Articulations

Anteriors		Posteriors			
Upper	Lower	A	B	C	D
M1	AL1	08	12	11	14
M2	AL2	08	13	11	14
M3	AL4	09	13	11	15
M4	AL7	09/10	13	11	15
R1	AL1	08	12	11	14
R2	AL2	07	13	11	14
R3	AL5	09	13	11	15
R4	AL6	09	13	11	15
R5	AL8	09	10	11	15
J1	AL1	08	12	11	14
J2	AL2	08	12	11	14
J3	AL3	09	13	11	15
J4	AL7	09/10	13	11	15
J5	AL8	09	10	11	15

Shades

Available in 16 shades and tones:

A1, A2, A3, A3.5, A4

B1, B2, B3, B4

C1, C2, C3, C4

D2, D3, D4

The dimensions are in millimeters.

Vertex Quint Teeth meet all the requirements of EN-ISO 22112:2006 and are not fluorescent.

Vertex Dental

Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist
 Postbus 10 | 3700 AA Zeist | The Netherlands
 Telefoon 030 - 69 76 749 | Fax 030 - 69 55 188
 info@vertex-dental.com | www.vertex-dental.com