

Vertex Dental

Casting System



Vertex™ Castagel Clear & Blue



Vertex™ Castagel can be used in combination with Vertex™ Castavaria or Vertex™ Castapress to produce effective a full or partial denture. With this system you use gel instead of plaster. Shrinkage will be prevented and time will be saved while transferring a wax setup by using a pouring type of acrylic.

Vertex™ Gelknife



With the Vertex™ Gelknife it is much easier to cut the gel into small pieces.

Vertex™ Castafask



Using the Vertex™ Castafask in combination with Vertex™ Castagel clear it is easy to control the whole process due to the transparency of the flask and the gel. The size is designed to control the shrinkage proces.

Vertex™ Canal cutters



The Vertex™ Canal cutter must be used in order to create pouring and ventilation channels in the gel.

Vertex™ Teeth holder



The Vertex™ Teeth holder is designed to release the wax from the teeth easily. Two types of teeth holders are available:

FVPOU0001 Vertex™ Teeth holder Metal



FVPOU0001A Vertex™ Teeth holder Plastic

The following products are related to the Casting System:



Vertex™ Thermometer

The covered thermometer can be mostly used for checking the Castagel's exact temperature before filling the flask.



Vertex™ PolyCure 25

For production due to it's capacity of 9 Vertex™ Castafasks with 25 litre chamber.



Vertex™ MultiCure

0° - 100°
0 - 6 bar
0 - 12 hrs
20 cm x 22 cm



Vertex™ GelMelt 8

Full-automatic gel melting machine with 8 litre tank.

Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist
P.O. box 10 | 3700 AA Zeist | The Netherlands
Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188
info@vertex-dental.com | www.vertex-dental.com