

Vertex Dental

High Gloss Assortment



Vertex-Dental made a selection of various high quality polishing brushes like; conventional polishing brushes, brushes with Scotch Brite, felt point, high gloss brushes and polishing paste.

Vertex™ Gloss brushes



A clean way of high gloss polishing without dust and loose wires coming out of the brushes during use, is possible with the long-life high polishing brushes from Vertex-Dental. Using these brushes with Vertex™ Gloss Paste will give excellent results. Vertex™ Brushes do not contain cotton, however microfibers and synthetic leather instead. Therefore improvements in term of quality, durability and user easiness are the user's benefits.

- FVPOL0011 Ø 100mm Polishing zeta microfiber 15 ply's
- FVPOL0013 Ø 100mm Felt cloth 6 ply's

Vertex™ L.C. Gloss varnish



Vertex™ L.C. Gloss Varnish is specially developed light curing varnish for the use on acrylics PMMA & MMA's. L.C Gloss varnish provides a high shining finish to acrylic areas that are hard to polish. Benefits: easy application, quick cure. Excellent bonding to hard acrylics.

Vertex™ High Gloss polishing paste



A clean way of polishing without dust and loose wires coming out of the brushes during use, is possible with the long-life high polishing brushes from Vertex-Dental. Using these brushes with Vertex™ Gloss Paste will give excellent results.

Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist
P.O. box 10 | 3700 AA Zeist | The Netherlands
Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188
info@vertex-dental.com | www.vertex-dental.com