

# Vertex Dental

## Machineries



### Vertex™ PolyCure 25



The Vertex™ PolyCure 25 has an adjustable water level, fixed pressure of 2,5 bar and a constant water temperature of 55 degrees Celsius. With this pressure vessel you are able to polymerize up to 9 Vertex™ Castaflasks in one run!

### Vertex™ MultiCure



- 20 cm x 22 cm = size of the chamber
- 0° - 100°
- 0 - 6 bar
- 0 - 12 hrs
- Water release valve

### Vertex™ GelMelt 8



This full-automatic gel melting machine has an 8 litre tank. By making use of microprocessor technology, very precise values can be realised. The melting and processing temperatures are programmed with push-buttons and potentiometers, and shown in the display.

### Vertex™ Eco Light Box



The Vertex™ Eco LightBox has a detachable tray with a black out front. A simple switch can be set at 5 or 10 minutes curing time. The high intensity and short polymerisation time, makes the machine suitable for polymerisation of individual light curing trays like the Vertex™ L.C. Trayplates.

The following products are related to the Machineries:



### Vertex™ Castaflask

Using in combination with Vertex™ Castagel clear it is easy to control the whole process due to the transparency of the flask and the gel.



### Vertex™ Metal flask

Upper, Lower and Clasp.



### Vertex™ Castagel Clear & Blue

For duplicate the waxmodel instead of plaster.



### Vertex™ LC Trayplates Pink, Blue & Transparent

Vertex™ Light Curing Trayplates with minimal residual layer after polymerisation makes processing easy. Each box contains 50 trays.

### Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist  
P.O. box 10 | 3700 AA Zeist | The Netherlands  
Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188  
info@vertex-dental.com | www.vertex-dental.com